



Bilirubin Meter ERMA B-111

Advanced, benchtop and precise instrument for measuring Total Bilirubin.



The advanced B111 Bilirubinmeter provides fast, accurate results for Total bilirubin level in neonates with only 40µl serum or plasma sample. It determines total bilirubin by dual wavelength photometric analysis. Suited for pediatric setups, small workload labs & hospitals because of its compact and portable design.

Clinical Implications

Bilirubin is a byproduct of the breakdown of hemoglobin by the liver. When hemoglobin is metabolized, insoluble bilirubin is produced. The insoluble (indirect) bilirubin is broken down into a water soluble (conjugated or direct) form. Direct bilirubin is then excreted. Most neonates have a degree of jaundice shortly after birth, because the liver is not fully developed. Newborn infants must be monitored to identify those who might develop severe hyperbilirubinemia.





ERMA B-111

Bilirubin Meter



Features :

Auto zero function:

- Manual blanking not required

Dual wavelength measurement system:

- 2 wavelengths 455nm & 575nm
- Plays major role to remove the interference caused due to Hemolysis & turbidity of sample

LED Light source:

- LED does not require replacement because it lasts long
- As compared to Lamp, LED are more suitable as there is no intensity fluctuations
- No warm up time required.

Flexible unit selection:

- Two units either mg/dl or umol/L are selectable, ease of reporting

Compact slot for Capillary:

- Accuracy of results are directly related to proper alignment of Capillary
- B 111 has easiest mode to ensure precise values.

Data transfer:

- External device can be connected to instrument through RS232 port
- USB Port



Specifications:

Light source	LED	Dimensions	150 (W) × 115 (H) × 190 (D) mm
Wavelength	455nm and 575nm	Net weight	1kg
Measuring range	0-30 mg/dl , 510 umol/L	Interface	RS232
Display	16x2 LCD module with backlight	Reading cuvette	Heparinized micro capillary holder
Auto calibration	Yes	Storage temperature	15-40°C
Sample type	Serum or plasma	Power supply	DC 9V(DC adapter)
Sample volume	40-80µL	Power consumption	<2W
Test time	<5sec		
Units	mg/dl or umol/L		

Accessories:

- External Adapter- 1pc
- Heparinized capillary- tube 1 box (100pcs)
- Sealing wax - 1pc
- Users manual - 1pc,

Heparinised micro capillary holder

Holds sample capillary in proper alignment which gives accurate results.

LCD Display

16 x 2 LCD module with backlight

Operation keys

Easy to operate



Di Distribusikan oleh:

PT TERAS SEKAWAN BERSAMA
JL.Radin Inten II No.7 Blok L-M, Rt.05/Rw.014
Kel.DurenSawit,Kec.DurenSawit
Jakarta Timur 13440

KEMENKES RI.NO.REG.AKL :.....

ERMA INC.Yushima 2-31-6,
Bunkyo-ku,Tokyo, Japan

